

Product Datasheet

Herpes Simplex virus, type 2 Antibody (biotin) (orb1273359)

Description	Herpes Simplex virus, type 2 Antibody (biotin)
Species/Host	Sheep
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, IF, WB
Immunogen	HSV type 2, strain G (human)
Preservatives	Herpes simplex virus, type 2 antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Herpes simplex virus, type 2 antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.