

Product Datasheet

Hepatitis C virus Antibody (orb1273375)

Description	Hepatitis C virus Antibody
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant (E. coli) polypeptide containing core + NS3 + NS4 regions
Preservatives	Hepatitis C virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Hepatitis C virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.