

Product Datasheet

Hepatitis B - surface antigen Antibody (biotin) (orb1273376)

Description	Hepatitis B - surface antigen Antibody (biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	WB
Immunogen	Subtypes ad & ay
Preservatives	Hepatitis B - surface antigen antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Hepatitis B - surface antigen antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.