

Product Datasheet

Hepatitis B - surface antigen Antibody (biotin) (orb1273380)

| | |
|--------------------------------|---|
| Description | Hepatitis B - surface antigen Antibody (biotin) |
| Species/Host | Goat |
| Reactivity | Virus |
| Conjugation | Biotin |
| Tested Applications | WB |
| Immunogen | Subtypes ad & ay |
| Preservatives | Hepatitis B - surface antigen antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added. |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Hepatitis B - surface antigen antibody can be stored at -20°C, stable for one year. |
| Note | For research use only |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |