

Product Datasheet

Hepatitis B - surface antigen Antibody (HRP) (orb1273381)

Description	Hepatitis B - surface antigen Antibody (HRP)
Species/Host	Goat
Reactivity	Virus
Conjugation	HRP
Tested Applications	WB
Immunogen	Subtypes ad & ay
Preservatives	Hepatitis B - surface antigen antibody is in a phosphate saline buffer with bovine serum albumin (10mg/ml) and preserved with 0.002% thimerosal.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Hepatitis B - surface antigen antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.