

## Product Datasheet

# Hepatitis B - surface antigen Antibody (FITC) (orb1273382)

<b>Description</b>	Hepatitis B - surface antigen Antibody (FITC)
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Virus
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Subtypes ad & ay
<b>Preservatives</b>	Hepatitis B - surface antigen antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Hepatitis B - surface antigen antibody can be stored at -20°C, stable for one year.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.