

Product Datasheet

Feline Leukemia Virus Antibody (orb1273403)

Description	Feline Leukemia Virus Antibody
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	purified FeLv p27
Preservatives	Feline Leukemia Virus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Feline Leukemia Virus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.