

Product Datasheet

Cytomegalovirus Antibody (FITC) (orb1273426)

Description	Cytomegalovirus Antibody (FITC)
Species/Host	Goat
Reactivity	Virus
Conjugation	FITC
Tested Applications	ICC, IF, WB
Immunogen	Purified virions of strain AD170
Preservatives	Cytomegalovirus antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Cytomegalovirus antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.