


Product Datasheet

Gabrg2 Antibody (orb1273537)

Description	Gabrg2 Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser327 of the GABAA receptor, $\gamma 2$ subunit.
Target	Gabrg2
Form/Appearance	Liquid
Concentration	batch dependent
Storage	GABAA Receptor (Ser327) antibody can be stored at -20°C and is stable at -20°C for at least 1 year.
Note	For research use only
Clonality	Polyclonal
MW	45
Uniprot ID	P18508
NCBI	P18508
Expiration Date	12 months from date of receipt.



Western blot of rat cortex showing speci...