
Product Datasheet

PLK1 Antibody (orb1273555)

Description	PLK1 Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Frog, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Thr210 of PLK1.
Target	PLK1
Form/Appearance	Liquid
Concentration	batch dependent
Storage	PLK1 antibody can be stored at -20°C and is stable at -20°C for at least 1 year.
Note	For research use only
Application notes	The antibody has been directly tested for reactivity in Western blots in rat and human tissues. It is anticipated that the antibody will react with bovine, canine, mouse, non-human primate, Xenopus and Zebrafish tissues based on the fact that these species have 100% homology with the amino acid sequence used as antigen.
Clonality	Polyclonal
MW	66
Uniprot ID	P53350
NCBI	P53350
Dilution Range	The antibody has been directly tested for reactivity in Western blots in rat and human tissues. It is anticipated that the antibody will react with bovine, canine, mouse, non-human primate, Xenopus and Zebrafish tissues based on the fact that these species have 100% homology with the amino acid sequence used as antigen.
Expiration Date	12 months from date of receipt.

