

---

## Product Datasheet

### Rph3a Antibody (orb1273567)

Description	Rph3a Antibody	
<b>Species/Host</b>	Rabbit	
<b>Reactivity</b>	Mouse, Rat	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	IHC, WB	Western blot of rat brain lysate showing...
<b>Immunogen</b>	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser234 of rabphilin 3A.	
<b>Target</b>	Rph3a	
<b>Form/Appearance</b>	Liquid	
<b>Concentration</b>	batch dependent	
<b>Storage</b>	Rabphilin antibody can be stored at -20°C and is stable at -20°C for at least 1 year.	
<b>Note</b>	For research use only	
<b>Application notes</b>	The antibody has been directly tested for reactivity in Western blots with rat tissue. It is anticipated that the antibody will react with mouse based on the fact that this species has 100% homology with the amino acid sequence used as antigen.	
<b>Clonality</b>	Polyclonal	
<b>MW</b>	82	
<b>Uniprot ID</b>	<a href="#">P47709</a>	
<b>NCBI</b>	<a href="#">P47709</a>	
<b>Dilution Range</b>	The antibody has been directly tested for reactivity in Western blots with rat tissue. It is anticipated that the antibody will react with mouse based on the fact that this species has 100% homology with the amino acid sequence used as antigen.	
<b>Expiration Date</b>	12 months from date of receipt.	

