


Product Datasheet

Kcnc1 Antibody (orb1273568)

Description	Kcnc1 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser503 of the voltage-gated potassium channel Kv3.1, conjugated to keyhole limpet hemocyanin (KLH).
Target	Kcnc1
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Potassium Channel antibody can be stored at -20°C and is stable at -20°C for at least 1 year.
Note	For research use only
Application notes	The antibody has been directly tested for reactivity in Western blots with rat and mouse tissue.
Clonality	Polyclonal
MW	100
Uniprot ID	P25122
NCBI	P25122
Dilution Range	The antibody has been directly tested for reactivity in Western blots with rat and mouse tissue.
Expiration Date	12 months from date of receipt.



IHC staining of medial nucleus of the tr...