


Product Datasheet

Synpr Antibody (orb1273599)

Description	Synpr Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to amino acid residues from the C-terminal region of rat synaptophysin 2, conjugated to keyhole limpet hemocyanin (KLH).
Target	Synpr
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Synaptophysin antibody can be stored at -20°C and is stable at -20°C for at least 1 year.
Note	For research use only
Clonality	Polyclonal
MW	38
Uniprot ID	P22831
NCBI	P22831
Expiration Date	12 months from date of receipt.



Western blot of rat hippocampal lysate showing specific binding of the anti-synaptophysin 2 protein.

Western blot of rat hippocampal lysate showing specific binding of the anti-synaptophysin 2 protein.