


Product Datasheet

Gabrb3 Antibody (orb1273652)

Description	Gabrb3 Antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phosphopeptide corresponding to amino acid residues surrounding the di-phospho-Ser408/409 of the rat GABAA receptor, $\beta 3$ subunit.
Target	Gabrb3
Form/Appearance	Liquid
Concentration	batch dependent
Note	For research use only
Clonality	Polyclonal
MW	53
Uniprot ID	P63079
NCBI	P63079
Expiration Date	12 months from date of receipt.



Western blot of rat hippocampal lysate s...