


Product Datasheet

Aqp2 Antibody (orb1273653)

| | |
|----------------------------|--|
| Description | Aqp2 Antibody |
| Species/Host | Rabbit |
| Reactivity | Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser269 of rat aquaporin 2. |
| Target | Aqp2 |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | 29/37 |
| Uniprot ID | P34080 |
| NCBI | P34080 |
| Expiration Date | 12 months from date of receipt. |



Western blot of vasopressin (AVP) treated rat kidney lysate showing specific immunolabeling of the ~29k and 37k glycosylated form of the AQP2 protein phosphorylated at Ser269.