


Product Datasheet

Aqp2 Antibody (orb1273654)

Description	Aqp2 Antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser264 of rat aquaporin 2.
Target	Aqp2
Form/Appearance	Liquid
Concentration	batch dependent
Note	For research use only
Clonality	Polyclonal
MW	29/37
Uniprot ID	P34080
NCBI	P34080
Expiration Date	12 months from date of receipt.

Anti-Phospho-Ser²⁶⁴ Aquaporin 2



Western blot of vasopressin (AVP) treated rat kidney lysate showing specific immunolabeling of the ~29k and 37k glycosylated form of the AQP2 protein phosphorylated at Ser²⁶⁴.

Western blot of vasopressin (AVP) treated...