


## Product Datasheet

### Ppp1r1b Antibody (orb1273655)

<b>Description</b>	Ppp1r1b Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser137 of rat DARPP-32.
<b>Target</b>	Ppp1r1b
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	32
<b>Uniprot ID</b>	<a href="#">Q6J4I0</a>
<b>NCBI</b>	<a href="#">Q6J4I0</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western blot of rat caudate lysate showing specific immunoblotting of rat caudate lysate immunoblotting of the ~32 kDa DARPP-32 phosphorylated at Thr<sup>137</sup> (Control). In the second lane immunoblotting is blocked by preadsorption of the phospho-peptide used as antigen.