


## Product Datasheet

### MECP2 Antibody (orb1273659)

<b>Description</b>	MECP2 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser80 of MECP2 (Methyl-CpG Binding Protein 2) protein.
<b>Target</b>	MECP2
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	75
<b>Uniprot ID</b>	<a href="#">P51608</a>
<b>NCBI</b>	<a href="#">P51608</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western blot of neonatal rat brain showing specific immunolabeling of the ~75k MECP2 protein phosphorylated at Ser<sup>80</sup>.