



Product Datasheet

ORF62 Antibody (orb1273675)

Description	ORF62 Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Ni-NTA purified recombinant human KSHV ORF62 expressed in E. Coli strain BL21 (DE3).
Target	ORF62
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	to KSHV ORF62 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:1:500 - 1,000ELISA:Determining optimal working dilutions by titration test.
Isotype	IgM
Clonality	Monoclonal
Uniprot ID	F5H8Y5
NCBI	F5H8Y5
Dilution Range	Western Blot:1:500 - 1,000ELISA:Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using KSHV ORF62 m...



Immunocytochemistry staining of TPA-induced B...