



Product Datasheet

ORF45 Antibody (orb1273676)

Description	ORF45 Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Ni-NTA purified recombinant human KSHV ORF45 expressed in E. Coli strain BL21 (DE3).
Target	ORF45
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	to KSHV ORF45 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:1:500 - 1,000ELISA:Determining optimal working dilutions by titration test.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	F5HDE4
NCBI	F5HDE4
Dilution Range	Western Blot:1:500 - 1,000ELISA:Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using KSHV ORF45 m...



Immunocytochemistry staining of TPA-induced ...