



Product Datasheet

K8 Antibody (orb1273716)

Description	K8 Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Ni-NTA purified recombinant human KSHV K8 α expressed in E. Coli strain BL21 (DE3).
Target	K8
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	KSHV K8 α monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:Dilution 1:500 - 1,000ICC:Dilution 1:50 - 1:100Determining optimal working dilutions by titration test.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q2HR82
NCBI	Q2HR82
Dilution Range	Western Blot:Dilution 1:500 - 1,000ICC:Dilution 1:50 - 1:100Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using KSHV K8alpha...



Immunocytochemistry staining of TPA-induced ...