



Product Datasheet

CIB1 Antibody (orb1273755)

Description	CIB1 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Ni-NTA purified truncated recombinant CIB1 expressed in E. Coli strain BL21 (DE3).
Target	CIB1
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	CIB1 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000IHC(F):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Clonality	Monoclonal
Uniprot ID	Q99828
NCBI	NP_001264693.1
Dilution Range	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000IHC(F):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using CIB1 monoclo...



Immunohistochemical analysis of paraffin...