

Product Datasheet

Crp Antibody (orb1273860)

Description	Crp Antibody
Species/Host	Gallus
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Rat C-Reactive protein peptide.
Target	Crp
Form/Appearance	Liquid
Concentration	batch dependent
Storage	2-8°C for one week and -20°C as an undiluted liquid for up to two years.
Note	For research use only
Application notes	This antibody can be used at 1:2000-1:5000 for both ELSIA and Western Blot analysis. Customers need to optimize the dilution when using the Abs in IHC or ICC.
Clonality	Polyclonal
Uniprot ID	P48199
NCBI	P48199
Dilution Range	This antibody can be used at 1:2000-1:5000 for both ELSIA and Western Blot analysis. Customers need to optimize the dilution when using the Abs in IHC or ICC.
Expiration Date	12 months from date of receipt.

