


Product Datasheet

MAPK8/MAPK9 Antibody (orb1273882)

Description	MAPK8/MAPK9 Antibody
Species/Host	Gallus
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide containing phosphorylated motif (pT PpY) of JNK
Target	MAPK8 / MAPK9
Preservatives	PBS
Form/Appearance	Liquid
Concentration	batch dependent
Storage	2-8°C for one week, -20°C as an undiluted liquid for up to two year.
Note	For research use only
Application notes	This antibody can be used at 1:1,000 to 1:4,000 for immunoblot testing and ELISA or 1:500-1:2000 for IHC or ICC.
Clonality	Polyclonal
Uniprot ID	P45983 , P45984
NCBI	P45984 , P45983
Dilution Range	This antibody can be used at 1:1,000 to 1:4,000 for immunoblot testing and ELISA or 1:500-1:2000 for IHC or ICC.
Expiration Date	12 months from date of receipt.



Detection of phosphorylated JNK/SAPK in ...