

## Product Datasheet

# SARS-CoV-2 (COVID-19) Spike ECD Monoclonal Antibody [B005] (HRP) (orb1273994)

<b>Description</b>	SARS-CoV-2 (COVID-19) Spike ECD Monoclonal Antibody [B005] (HRP)
<b>Species/Host</b>	Human
<b>Reactivity</b>	Virus
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, NeA, WB
<b>Immunogen</b>	S-ECD recombinant protein
<b>Target</b>	S
<b>Preservatives</b>	PBS with 0.02% sodium azide, 10% glycerol, pH 7.2
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: 1:1000~2000
<b>Clonality</b>	Monoclonal
<b>NCBI</b>	<a href="#">QHD43416.1</a>
<b>Dilution Range</b>	WB: 1:1000~2000
<b>Expiration Date</b>	12 months from date of receipt.