


Product Datasheet

SARS-CoV-2 (COVID-19) Spike RBD Monoclonal Antibody [B001] (Biotin) (orb1274011)

Description	SARS-CoV-2 (COVID-19) Spike RBD Monoclonal Antibody [B001] (Biotin)
Species/Host	Human
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, NeA, WB
Immunogen	S-RBD recombinant protein
Target	S
Preservatives	PBS with 0.02% sodium azide, 10% glycerol, pH 7.2
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	WB: 1:1000~2000
Clonality	Monoclonal
NCBI	QHD43416.1
Dilution Range	WB: 1:1000~2000
Expiration Date	12 months from date of receipt.



Western Blot Validation with Recombinant...