


Product Datasheet

SARS-CoV-2 (COVID-19) Spike RBD Monoclonal Antibody [B001] (azide free) (orb1274012)

Description	SARS-CoV-2 (COVID-19) Spike RBD Monoclonal Antibody [B001] (azide free)
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, NeA, WB
Immunogen	S-RBD recombinant protein
Target	S
Preservatives	PBS, pH 7.2
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	ELISA: 1:320000; WB: 1:1000~2000
Clonality	Monoclonal
NCBI	QHD43416.1
Dilution Range	ELISA: 1:320000; WB: 1:1000~2000
Expiration Date	12 months from date of receipt.



Western Blot Validation with Recombinant...

ELISA experiment

	1	2	
A	2.347	2.399	1:1000
B	2.036	2.216	1:10000
C	2.297	2.139	1:20000
D	2.441	2.624	1:40000
E	2.104	2.262	1:80000
F	1.794	1.713	1:160000
G	1.186	1.213	1:320000
H	0.105	0.066	Blank control

ELISA Test.
Coating original concentration...