


## Product Datasheet

# SARS-CoV-2 (COVID-19) Nucleocapsid Monoclonal Antibody [N008] (orb1274021)

<b>Description</b>	SARS-CoV-2 (COVID-19) Nucleocapsid Monoclonal Antibody [N008]
<b>Species/Host</b>	Human
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Nucleocapsid recombinant protein
<b>Target</b>	N
<b>Preservatives</b>	PBS with 0.02% sodium azide, 10% glycerol, pH 7.2
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA: 1:40000; WB: 1:2000~3000
<b>Clonality</b>	Monoclonal
<b>NCBI</b>	<a href="#">QHD43423.2</a>
<b>Dilution Range</b>	ELISA: 1:40000; WB: 1:2000~3000
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot Validation with Recombinant...

ELISA experiment

	1	
A	2.171	1:10000
B	1.73	1:20000
C	1.271	1:40000
D	0.829	1:80000
E	0.497	1:160000
F	0.382	1:320000
G	0.28	1:640000
H	0.137	Blank control

ELISA Test. Coating original concentrati...