


## Product Datasheet

# SARS-CoV-2 (COVID-19) Nucleocapsid Antibody (NT) (orb1274051)

<b>Description</b>	SARS-CoV-2 (COVID-19) Nucleocapsid Antibody (NT)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide located within the N terminus of SARS-CoV-2 (COVID-19) Nucleocapsid.
<b>Target</b>	N
<b>Preservatives</b>	PBS with 0.02% sodium azide, 10% glycerol, pH 7.2
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: 1:500-1000
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>NCBI</b>	<a href="#">QHD43423.2</a>
<b>Dilution Range</b>	WB: 1:500-1000
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot Validation with Recombinant...

ELISA experiment

	1	
A	2.133	1:100
B	2.103	1:200
C	1.445	1:500
D	1.037	1:1000
E	0.494	1:5000
F	0.52	1:10000
G	0.469	1:20000
H	0.167	Blank control

ELISA Test. Coating original concentration...