

Product Datasheet

Orf1ab Antibody (orb1274056)

| | |
|----------------------------|--|
| Description | Orf1ab Antibody |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Synthetic peptide located within NSP1 of SARS-CoV-2 (COVID-19) pp1ab. |
| Target | Orf1ab |
| Preservatives | PBS with 0.02% sodium azide, 10% glycerol, pH 7.2 |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Application notes | WB: 1:500-1000 |
| Isotype | IgG |
| Clonality | Polyclonal |
| NCBI | QHD43415.1 |
| Dilution Range | WB: 1:500-1000 |
| Expiration Date | 12 months from date of receipt. |