

## Product Datasheet

# MERS Coronavirus (HSZ-Cc) Membrane Antibody (orb1274121)

<b>Description</b>	MERS Coronavirus (HSZ-Cc) Membrane Antibody
<b>Species/Host</b>	Sheep
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	GST-MERS CoV HSZ-Cc M Protein
<b>Target</b>	M
<b>Preservatives</b>	Antibody supplied in PBS.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Aliquot and store at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">A0A2H4TED7</a>
<b>Expiration Date</b>	12 months from date of receipt.