

Product Datasheet

SARS Coronavirus (HSZ-Cc) Membrane Antibody (orb1274122)

Description	SARS Coronavirus (HSZ-Cc) Membrane Antibody
Species/Host	Sheep
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	GST-SARS CoV HSZ-Cc M Protein
Target	M
Preservatives	Antibody supplied in PBS.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Aliquot and store at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P59596
Expiration Date	12 months from date of receipt.