

Product Datasheet

MERS Coronavirus (HSZ-Cc) Spike Antibody (orb1274131)

Description	MERS Coronavirus (HSZ-Cc) Spike Antibody
Species/Host	Sheep
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	GST-MERS CoV HSZ-Cc Spike
Target	S
Preservatives	Antibody supplied in PBS.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Aliquot and store at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	W6A028
Expiration Date	12 months from date of receipt.