


Product Datasheet

SARS Nucleocapsid Antibody [14B3D] (orb1274141)

Description	SARS Nucleocapsid Antibody [14B3D]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS N protein (no tag)
Target	N
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	≥1 mg/mL
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	WB(1:500-2000), ELISA(1:5000-1000)
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P59595
NCBI	P59595
Dilution Range	WB(1:500-2000), ELISA(1:5000-1000)
Expiration Date	12 months from date of receipt.



Western Blot Validation of SARS Nucleoca...