


## Product Datasheet

### SARS Nucleocapsid Antibody [5B1B] (orb1274144)

<b>Description</b>	SARS Nucleocapsid Antibody [5B1B]
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant SARS N protein (no tag)
<b>Target</b>	N
<b>Preservatives</b>	0.01M PBS(pH7.4) with 0.1% Proclin300
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	≥1 mg/mL
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	WB(1:500-2000), ELISA(1:5000-1000)
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P59595</b>
<b>NCBI</b>	<b>P59595</b>
<b>Dilution Range</b>	WB(1:500-2000), ELISA(1:5000-1000)
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot Validation of SARS Nucleoca...