

## Product Datasheet

### orf1ab Antibody (HRP) (orb1274147)

<b>Description</b>	orf1ab Antibody (HRP)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC-P
<b>Immunogen</b>	KLH synthetic peptide derived from SARS ORF1a polyprotein
<b>Target</b>	orf1ab
<b>Preservatives</b>	1% rAlbumin, 50% glycerol and 0.09% Gentamicin
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA(1:500-1000), IHC-P(1:200-400)
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P0C6X7</a>
<b>NCBI</b>	<a href="#">P0C6X7</a>
<b>Dilution Range</b>	ELISA(1:500-1000), IHC-P(1:200-400)
<b>Expiration Date</b>	12 months from date of receipt.