

Product Datasheet

orf1ab Antibody (Cy3) (orb1274153)

Description	orf1ab Antibody (Cy3)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Cy3
Tested Applications	IF
Immunogen	KLH synthetic peptide derived from SARS ORF1a polyprotein
Target	orf1ab
Preservatives	1% rAlbumin, 50% glycerol and 0.09% sodium azide
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	IF(IHC-P)(1:50-200)
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C6X7
NCBI	P0C6X7
Dilution Range	IF(IHC-P)(1:50-200)
Expiration Date	12 months from date of receipt.