

Product Datasheet

orf1ab Antibody (orb1274154)

Description	orf1ab Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P
Immunogen	KLH synthetic peptide derived from SARS ORF1a polyprotein
Target	orf1ab
Preservatives	1% rAlbumin, 50% glycerol and 0.09% sodium azide
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	ELISA(1:500-1000), IHC-P(1:200-400), IF(IHC-P)(1:50-200)
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C6X7
NCBI	P0C6X7
Dilution Range	ELISA(1:500-1000), IHC-P(1:200-400), IF(IHC-P)(1:50-200)
Expiration Date	12 months from date of receipt.