

Product Datasheet

SARS Spike Glycoprotein S1 Antibody (orb1274155)

Description	SARS Spike Glycoprotein S1 Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH synthetic peptide derived from SARS spike glycoprotein S1 (101-200/1255)
Target	S
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	ELISA(1:500-1000)
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P59594
NCBI	P59594
Dilution Range	ELISA(1:500-1000)
Expiration Date	12 months from date of receipt.