

Product Datasheet

SARS Spike Antibody (Biotin) (orb1274157)

Description	SARS Spike Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH synthetic peptide derived from SARS S-protein
Target	S
Preservatives	1% rAlbumin, 50% glycerol and 0.09% sodium azide
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	WB(1:300-1000), ELISA(1:500-1000), IHC-P(1:200-400)
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P59594
NCBI	P59594
Dilution Range	WB(1:300-1000), ELISA(1:500-1000), IHC-P(1:200-400)
Expiration Date	12 months from date of receipt.