

Product Datasheet

SARS Spike Antibody (Cy5.5) (orb1274160)

Description	SARS Spike Antibody (Cy5.5)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Cy5.5
Tested Applications	IF
Immunogen	KLH synthetic peptide derived from SARS S-protein
Target	S
Preservatives	1% rAlbumin, 50% glycerol and 0.09% sodium azide
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	IF(IHC-P)(1:50-200)
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P59594
NCBI	P59594
Dilution Range	IF(IHC-P)(1:50-200)
Expiration Date	12 months from date of receipt.