

## Product Datasheet

### SARS Spike Antibody (Cy5) (orb1274161)

<b>Description</b>	SARS Spike Antibody (Cy5)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Cy5
<b>Tested Applications</b>	IF
<b>Immunogen</b>	KLH synthetic peptide derived from SARS S-protein
<b>Target</b>	S
<b>Preservatives</b>	1% rAlbumin, 50% glycerol and 0.09% sodium azide
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	IF(IHC-P)(1:50-200)
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P59594</a>
<b>NCBI</b>	<a href="#">P59594</a>
<b>Dilution Range</b>	IF(IHC-P)(1:50-200)
<b>Expiration Date</b>	12 months from date of receipt.