

Product Datasheet

ACE2 Antibody (orb1274170)

Description	ACE2 Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	C terminal human Angiotensin Converting Enzyme 2
Target	ACE2
Preservatives	1xTBS (pH7.4), 1%rAlbumin, 40%Glycerol and 0.05% Sodium Azide
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	WB(1:300-1000), IF(ICC)(1:50-200), IHC()
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q9BYF1
NCBI	Q9BYF1
Dilution Range	WB(1:300-1000), IF(ICC)(1:50-200), IHC()
Expiration Date	12 months from date of receipt.