


---

## Product Datasheet

# SARS-CoV-2 (COVID-19) Spike RBD Single Domain Antibody [T4P3-B7] (orb1274182)

**Description**

SARS-CoV-2 (COVID-19) Spike RBD Single


**Species/Host**

Camelus

**Reactivity**

Virus

**Conjugation**

Unconjugated

**Tested Applications**

ELISA

 ELISA  
Validation  
with  
SARS-  
CoV-2  
(COVID-...

**Immunogen**

SARS-CoV-2 S protein RBD containing C-terminal His Tag. The protein was expressed in human 293 cells (HEK293). It contains amino acids Arg 319 - Lys 537.

**Target**

S

**Preservatives**

SARS-CoV-2 (COVID-19) S1 RBD Antibody [T4P3-B7] Antibody is supplied in PBS.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

SARS-CoV-2 (COVID-19) S1 RBD antibody should be stored in working aliquots at -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Isotype**

sdAb

**Clonality**

Recombinant

**NCBI**
[QHD43417](#)
**Expiration Date**

12 months from date of receipt.