



Product Datasheet

PDCD1 Antibody (orb1274202)

Description	PDCD1 Antibody
Species/Host	Camelus
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	PD-1 antibody was raised against a recombinant protein corresponding to amino acids 25 - 167 of human PD-1.
Target	PDCD1
Preservatives	PD-1 Antibody is supplied in PBS.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	PD-1 antibody should be stored in working aliquots at -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	PD-1 antibody can be used for detection of PD-1 by ELISA and FACS (10 µg/ml).
Isotype	sdAb
Clonality	Recombinant
Uniprot ID	Q15116
NCBI	NP_005009
Dilution Range	PD-1 antibody can be used for detection of PD-1 by ELISA and FACS (10 µg/ml).
Expiration Date	12 months from date of receipt.



Titration
ELISA of PD-1.



Flow
cytometry
analysis of
PD-1
transfec...