



Product Datasheet

YBX1 Antibody (orb1274422)


Description	YBX1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human YBX1 (NP_004550.2).
Target	YBX1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IHC: 1:100 - 1:200IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 49kDa
Uniprot ID	P67809
Dilution Range	WB: 1:500 - 1:2000IHC: 1:100 - 1:200IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Hep...



Western blot analysis of extracts from n...



Immunohistochemistry of paraffin-embedde...