



Product Datasheet

GSK3B Antibody (orb1274437)


Description	GSK3B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human GSK3B (NP_001139628.1).
Target	GSK3B
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000IF: 1:50 - 1:100IP: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 46kDa
Uniprot ID	P49841
Dilution Range	WB: 1:500 - 1:1000IF: 1:50 - 1:100IP: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Western blot analysis of extracts from n...



Immunofluorescence analysis of HeLa cell...