

Product Datasheet

Recombinant SARS-CoV-2 (COVID-19) S2 antibody [AbD103] (orb1274577)

| | |
|----------------------------|--|
| Description | Recombinant SARS-CoV-2 (COVID-19) S2 antibody [AbD103] |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant SARS-CoV-2 S2 protein |
| Preservatives | PBS, pH7.5 |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage. |
| Note | For research use only |
| Application notes | WB:1:1000-2000ELISA:1:5000-1:10000 |
| Isotype | Human IgG |
| Clonality | Recombinant |
| Dilution Range | WB:1:1000-2000ELISA:1:5000-1:10000 |
| Expiration Date | 12 months from date of receipt. |