


## Product Datasheet

# Recombinant CoV-2/COVID-19 Spike Glycoprotein Antibody (orb1274596)

|                            |                                                                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Recombinant CoV-2/COVID-19 Spike Glycoprotein Antibody                                                                                                                      |
| <b>Species/Host</b>        | Human                                                                                                                                                                       |
| <b>Reactivity</b>          | Virus                                                                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                       |
| <b>Immunogen</b>           | Recombinant SARS-CoV-2/COVID-19 Spike Glycoprotein                                                                                                                          |
| <b>Preservatives</b>       | PBS,pH7.5                                                                                                                                                                   |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                      |
| <b>Concentration</b>       | batch dependent                                                                                                                                                             |
| <b>Storage</b>             | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage. |
| <b>Note</b>                | For research use only                                                                                                                                                       |
| <b>Application notes</b>   | ELISA:1:5000-1:10000                                                                                                                                                        |
| <b>Isotype</b>             | Human IgM                                                                                                                                                                   |
| <b>Clonality</b>           | Recombinant                                                                                                                                                                 |
| <b>Dilution Range</b>      | ELISA:1:5000-1:10000                                                                                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                             |



Coomassie blue staining non-reduced and ...