



Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [85B4] (orb1274603)

Description	SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [85B4]
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	SARS-CoV-2 nucleoprotein (NP) specific peptide
Target	N
Preservatives	1X PBS, 0.02% Na3
Form/Appearance	Liquid
Concentration	batch dependent
Storage	SARS-CoV-2 (COVID-19) Nucleoprotein Antibody [85B4] can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	ELISA
Isotype	IgG1k
Clonality	Recombinant
NCBI	QHD43423
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Microtiter wells were coated with SARS-CoV-2 Nucleocapsid Protein (1 µg/ml).



Western blot analysis of SARS-CoV2 and SARS-CoV.